
Etac Nova Stand Up

molift[®]
by Etac


User Manual - Bruksanvisning - Brukermanual - Brugervejledning - Käyttö-
pas - Bedienungsanleitung - Gebruikershandleiding - Manuel d'utilisation
- Manuale per l'utente - Manual de usuario


BM28199 Rev. B 2021-02-09

MD

CE

Contents

User manual - English - 4
Bruksanvisning - Svenska - 8
Brukermanual - Norsk - 12
Brugervejledning - Dansk - 16
Käyttöohje - Suomi - 20
Gebrauchsanweisung - Deutsch - 24
Handleiding - Nederlands - 28
Manual d'utilisation - Français - 32
Manual utente - Italiano - 36
Manual de usuario - Español - 40

General

Declaration of conformity

The Etac Nova StandUp and related accessories described in this operator manual are CE marked in accordance with EU Council Directive MDR (EU) 2017/745 concerning medical devices, class 1, and has been tested and approved by a third party according to standard EN ISO 10535:2006.

Conditions for Use: Lift and transfer of a person will always pose a certain risk and only trained personnel should use the equipment and accessories covered by this user manual.

Warranty: 2-year warranty against defects in workmanship and materials of our products. For Terms and conditions, visit www.etac.com

Visit www.etac.com for download of documentation to ensure you have the latest version.

Important

This User Manual contains important safety instructions and information regarding the use of the sling and accessories. In this manual the user is the person being lifted. The assistant is the person operating the hoist.

This symbol indicates important information related to safety. Follow these instructions carefully.

Read User Manual before use!
It is important to fully understand the content of the user manual before attempting to use the equipment. Read the manual for both hoist and sling

General Safety Precautions: Only use accessories and slings that are adjusted to fit the user, type of disability, size, weight and type of transfer. It is important that the sling has been tested with the individual user and for the intended lifting situation. Make a decision on whether one or more assistants are required. Plan the lifting operation in advance to ensure that it is as safe and smooth as possible. Remember to work ergonomically.

Assess the risk and take notes.

You as the assistant are responsible for the safety of the user!

Molift slings should only be used to lift people. Never use the sling to lift or move objects of any kind.


If maximum load (SWL) differs between hoist, suspension and body support unit, then the lowest maximum load shall always be used

The sling should not be stored in direct sunlight.

Labels and symbols

1. Main label
2. Size label
3. Batch and serial number label w. GS1-128 barcode.
4. Product label
5. Periodic inspection label/ Name tag

Symbols:

Manufacturer

Date of manufacture.
YYYY-MM-DD
(year/month/date)

CE marked

Refer to user manual

This way up, this side is outside

Max user weight

Medical Device

Wash symbols:

Slings can be washed at temperatures ranging from 60 up to 85° C

Tumble drying, max 60° C.

Do not iron


Do not dry clean


Do not bleach

About Etac Nova StandUp

Etac supplies a wide selection of slings for different types of transfers. The Nova StandUp is designed to distribute weight and pressure evenly to make the transfer as comfortable as possible. The Etac Nova Stand Up is intended to be used with the sit-to-stand hoist Etac Nova 500 NG. The Etac Nova StandUp is available in sizes XS – XXL, in padded polyester. The user can move from one sitting position to another, the product is suitable for transfer to a toilet. Etac Nova Standup is intended for users who can follow instructions, have good stability in their upper body and can weight bear on at least one leg, but need assistance when standing. Nova StandUp has a curved cut back, is padded to give good comfort and requires only small individual adjustments. It is placed around the body to provide support to the lower back. The sling has an anti-slip lining inside and a waist belt, which helps keep the sling in place during transfers. The sling has a lifting strap on each side which attaches to the suspension on the hoist. See the combination list for the correct sling and suspension.

Technical data

Safe Working Load (SWL):

XS-XXL: 200 kg (441 lbs)

Material:

Polyester, Polyethylene



Size	A	B
XS	30	69
S	35	82
M	35	93
L	35	119
XL	35	141
XXL	35	158

Expected Lifetime:

The product has an expected service life of 1 to 5 years under normal use. The service life of the product varies depending on usage frequency, materials, loads and how often it is washed. Washing at higher temperatures wears out the material faster.

Combination list

Etac Nova StandUp is designed for use with sit-to-stand hoist Etac Nova 500 NG.

	Sling size		
	XS-S	M-L	XL-XXL
 Etac Nova 500 NG	✓	✓	✓
 Etac Nova 500 NG with extra bar	✓	✓	✓

Before use / Daily check

Before using a sling for the first time it should be marked with the date of first time use on the periodic inspection label.

Inspection to be performed daily or before use:

- Make sure the user has the required ability for the sling
- Make sure the sling suitable for the hoist/ suspension to be used.
- Make sure the sling does not have visible damage or frays.

Periodic inspection

Periodic inspection of the non-rigid body support unit should be undertaken at the time intervals stated by the manufacturer, but at least every 6 months. More frequent inspections may be required where a non-rigid body support unit is used or cleaned more frequently than normal


Never use a faulty or damaged sling as it can break and cause personal injury. Destroy and discard damaged and old slings

Make sure not to damage or remove labels when cleaning.

Size guide


Measurements in the table should be used as a reference. Correct size depends on the user's weight, function and body shape. Measure the waist size and back height of the user in sitting position. Try the sling with most the correct measurements to make sure it fits.

Size	User weight - kg (lbs)	Waist (W)
XS	17-25 (37-55)	60-75cm
S	25-50 (55-110)	70-95cm
M	45-95 (99-210)	90-105cm
L	90-165 (198-353)	100-115cm
XL	160-240 (353-529)	110-130cm
XXL	230-300 (507-660)	125-145cm


Sitting position - Sling application


1. Ask the user to lean forward, support them if necessary with one hand. Slide the sling down behind the user's back. Place the sling approximately at waist height on the back.

2. Place the belt around the waist of the user, connect the buckle and tighten the belt.


3. Place the hoist in position for fastening the sling to the suspension. Place the user's feet on the footplate and fasten the strap on the leg support around the user's leg.


4. Fasten the sling to the suspension. Make sure the lifting straps are fastened with equal length on both sides.

Transfer

1. Check that the sling is correctly fitted around the user and that the strap loops are correctly fitted to the suspension hooks.
2. Make sure the user has both feet placed firmly on the footplate on the hoist and the strap on the leg support is properly fastened around the user's leg.
3. Start lifting until the lifting straps are stretched without lifting the user. Ensure that both loops of the sling are securely fastened to avoid the user slipping or falling out of the sling.
4. Make sure the buckle is properly fastened
5. Lift the user to standing position. Tighten and adjust the waist belt if necessary.

The user must lean back slightly when being lifted, and if possible hold onto the handles. The waist belt should be retightened when lifting begins and when the sling is loaded.

Never lift a user higher than is necessary to carry out a sit to stand transfer. Remember that wheels on a mobile hoist must NOT be locked during the transfer

Be careful during movement, the user may swing during turns, stops and starts. Be careful when manoeuvring close to furniture and similar to ensure the suspended user does not hit anything.

Avoid deep pile carpets, high thresholds, uneven surfaces or other obstacles that may block the castors. The hoist may become unstable if forced over such obstacles increasing the risk of it tipping over.


Never leave a user unattended while being hoisted.

The hoist shall not be used to lift or move users on sloping surfaces

Sitting position - Removing the sling

1. Lower the user into a sitting position. Release the lifting straps from the suspension. Remove the hoist, and unlock the waist belt buckle.

2. Stand next to the chair and ask the user to lean forward. Support the user with one hand if necessary. Pull the sling away with one hand, and ask the user to lean back in the chair.

Remove the sling gently to ensure the user does not fall forward.

Recycling

Recycling instructions are available at Etac.com

Troubleshooting

Symptom	Possible Cause/Action
The sling does not fit the user.	Wrong sling size. Try other size.
	The sling is not correctly positioned. Reposition the sling.
The sling is sliding up to the users armpits	The user is not leaning backwards. Make the user lean backwards when lifting.
	The waist belt is loose. Tighten the waist belt properly.
	The user is not able to take load on their legs/feet. Evaluate whether a sit to stand hoist is right or if the user requires a passive hoist.
User is sliding through/out of the sling.	The lifting straps are not fastened correctly. Make sure both lifting straps are equal in length.
	The sling is too big. Try a different size
	The user is not leaning backwards. Make the user lean backwards when lifting.
	The waist belt is loose. Tighten the waist belt properly.
	The user is wearing slippery clothes. Change or remove some of the clothing.
	The user is not able to weight bear on their legs/feet. Evaluate whether a sit to stand hoist is right or if the user requires a passive hoist.

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Etac A/S