Hill-Rom. Liko_m MultiStrap_m





For a Wide Range of Use in Healthcare.

MultiStrap^{••} is a lifting accessory from Liko^{••} which completes the other lifting and repositioning aids in Liko's wide range. MultiStrap^{••} has been developed to meet a variety of needs in care and can be used to hold arms, legs and other parts of the body, for example, during treatment or washing.

Together with a patient lift, MultiStrap[™] facilitates many heavy lifting tasks such as wound care (e.g. burn injuries), hygiene procedures, treatment (e.g. catheterizing or dressing wounds), examinations, placement of bedpans and dressing/undressing. The MultiStraps can also be used to turn a patient to the side. MultiStrap[™] can be used individually or in pairs. For turning in bed, we recommend using it in pairs with one of the wider models.

Important Features

- Facilitates lifting/holding e.g. arms and legs during treatment/washing
- Three washable MultiStraps in different sizes
- Also available are MultiStraps in disposable fabric intended for short-term use by an individual patient



Enhancing Outcomes for Patients and Their Caregivers...

Different Models and Areas of Use.

Liko MultiStrap Standard

Width 10 cm (4 inches) Liko₁₁ MultiStrap₁₁ Standard is usually the easiest to apply and allows better access to a larger area of skin. It is especially suitable for lifting and holding an arm or leg.

Liko MultiStrap Wide

Width 25 cm (10 inches) Liko₁₁ MultiStrap₁₁ Wide has a greater support area and distributes pressure better thanks to its reinforcement. It can be used to lift a patient's leg or behind, or to turn a patient in bed.

Liko MultiStrap Ultra

Width 40 cm (16 inches) Developed especially for bariatric patients but can also be used to turn lighter patients in bed. Has reinforcement to distribute pressure more evenly.



A FixStrap 135 cm (54 in) is included with all Liko^{III} MultiStraps and is used to fix one end of the MultiStrap. to the bed when turning a patient.

Read the relevant instruction guides for further information on areas of use and practical application. You can also contact your closest Likom representative for information and advice.

Product Information

Product Designation Liko "MultiStrap "Standard, incl. FixStrap Liko MultiStrap Wide, incl. FixStrap Liko "MultiStrap "Ultra, incl. FixStrap

Product No. Maximum Load 200 kg (440 lbs) 3695115 200 kg (440 lbs) 3695116 3695117 200 kg (440 lbs) Product No. 3606001

Accessories

FixStrap 135 cm (54 in) Spatula 600 IC

3684002 (can facilitate placement of Liko "MultiStrap. Standard and Wide) Fabric Polyester Polyester 25 cm (10 in) Polyester 40 cm (16 in)

Width 10 cm (4 in) MultiStrap wis a class 1 medical device product.

Design and Quality by Liko in Sweden.

Liko is quality-certified in accordance with ISO 9001 and its equivalent for the medical device industry, ISO 13485. Liko is certified in accordance with environmental standard ISO 14001. C E 间

Hill-Rom is a trademark or registered trademark of Hill-Rom Services Inc. Liko is a trademark or registered trademark of Liko R&D AB. MultiStrap is a trademark or registered trademark of Liko R&D AB.



Enhancing Outcomes for Patients and Their Caregivers....

+33 (0)2 97 50 92 12
+44 (0)1530 411000
+49 (0)211 16450-0
+31 (0)347 32 35 32
+ 39 (0)2 950541
+41 (0)21 706 21 30
+41 (0)21 706 21 38

Österreich	+43 (0)2243 28550
Ireland	+353 (0)1 413 6005
Iberia	+34 (0)93 6856000
Nordic Countries	+46 (0)20 78 10 30
distributor or conta	ontact your local Hill-Rom ct your Area Manager via +1 812 934 8173

www.hill-rom.com www.liko.com

Hill-Rom_® reserves the right to make changes without notice in design, specifications and models. The only warranty Hill-Rom_® makes is the express written warranty extended on the sale or rental of its products. © 2012 Hill-Rom_® SARL, ALL RIGHTS RESERVED.