


Fitting no.1

Short and straight fitting. Well suited if the chair is very slim and the side supports are supposed to be close to the body


Measurements
H:6 x L:10½ cm



Fitting no.2

Short and straight fitting. Well suited if the chair is very slim and the side supports are supposed to be close to the body


Measurements
H:6 x L:13½ cm



Fitting no.3

Straight fitting in stainless steel. Used on bathing chairs


Measurements
H:4½ x L:12 cm



Fitting no.4

Small fitting with an angle on the plate, so that it fits for the curved back of the chairs


Measurements
H:5½ x L:12 cm x 7½°



Fitting no.5

Stainless steel fitting, curved. The rotation axis is displaced backward. Protects children against hurting their arms


Measurements
H:5½ x L:13 cm



Fitting no.6

Curved fitting with the rotation axis displaced backward. Protects children against hurting their arms


Measurements
H:4½ x L:12 cm x 7½°



Fitting no.7

With this fitting the side support cannot be adjusted in width, which makes it important to choose the correct bar


Measurements
H:8 x L:7 cm



Fitting no.8

With this fitting the side support can not be adjusted in width, which makes it important to choose the correct bar


Measurements
H:8 x L:7 cm



Fitting no.9

Curved fitting. The rotation axis is displaced backward. Protects the child against hurting its arms


Measurements
H:4½ x L:12 cm



Curved cushion no.1

The curved cushions are used when sideways and forward support is needed


Size **Measurements**
one H:9 x W:19 cm



Curved cushion no.2

The curved cushions are used when sideways and forward support is needed


Size **Measurements**
one H:6½ x W:17 cm



Curved cushion no.3

The curved cushions are used when sideways and forward support is needed


Size **Measurements**
one H:6½ x W:15 cm



Standard cushion no.4

The standard cushions are only used for sideways support. They provide the user with good freedom of movement of the arms


Size	Measurements
one	H:15½ x W:12 cm



Standard cushion no.5

The standard cushions are only used for sideways support. They provide the user with good freedom of movement of the arms


Size	Measurements
one	H:10 x W:10 cm



Standard cushion no.6

The standard cushions are only used for sideways support. They provide the user with good freedom of movement of the arms


Size	Measurements
one	H:5½ x W:11½ cm



Bar, straight

This bar is used to gain the largest width between the cushions


Size	Measurements
no.10	L:12 cm
no.11	L:16 cm
no.15	L:9½ cm
no.16	L:16 cm



Bar, curved, 25 mm

This bar is used when the back of the user is a little smaller than the back of the chair


Size	Measurements
no.25	L1:12 x L2:2½ cm
no.26	L1:16 x L2:2½ cm



Bar, curved, 45 mm

This bar is used when the back of the user is a little smaller than the back of the chair

Size	Measurements
no.45	L1:12 x L2:4½ cm
no.46	L1:16 x L2:4½ cm



Bar, curved, 65 mm

This bar is used to gain the smallest width between the cushions

Size	Measurements
no.65	L1:12 x L2:6½ cm
no.66	L1:16 x L2:6½ cm