



Molift UnoSling

Prevents transmission of infection





Molift UnoSling HighBack

A patient specific sling for use where there are strict standards for infection control.

High requirements for infection control

The Molift UnoSling HighBack is a personal disposable product designed to meet the demands and challenges of today's hospital environment. Slings for short-term use are specific to a single user, and an effective solution for preventing the transmission of infection.

The UnoSling has a high back and is suitable for most users in the most common hoisting situations, such as between bed and wheelchair or to the toilet. The sling remains with the user during the entire length of their stay in hospital.

Sling for short-term use

The UnoSling is manufactured in a strong non-woven material which is not washable. The sling is discarded if it becomes dirty or damaged, or if the user no longer needs it.



Molift UnoSling HighBack	Number in packaging	User weight	Art. no.	Material
S	20	25-50 kg	3060020	Polypropylene and polyester. Disposable product Do not wash. The UnoSling should be discarded when the user no longer needs it, or if it becomes damaged or dirty.
M	20	45-95 kg	3060030	
L	20	90-160 kg	3060040	
XL	20	160-240 kg	3060050	
XXL	10	230-300 kg	3060060	



max 300 kg
max 660 lbs



Details



Disposable product

The Molift UnoSling is supplied in packaging which makes it easy to store and take out a new sling when needed.



Suitable for use with 2-point or 4-point slingbar

The Molift UnoSling is suitable for use with either 2-point or 4-point suspension.



Label with user's name

Large, clear label with practical information. The size is shown clearly and it is easy to add the user's name, which ensures safe use.



Not washable

Do not wash. The UnoSling should be discarded when the user no longer needs it, or if it becomes damaged or dirty.





“All people should have the opportunity to create their own lives and their own daily living on their own terms, irrespective of their physical conditions.”

You can also find news and updated information at www.etac.se.

Etac AB
Kista Science Tower
SE-164 51 Kista, Sweden
Tel +46 371 - 58 73 30
Fax +46 371 - 58 73 90
info@etac.se www.etac.com



Mediq Suomi Oy
PL 115, 02201 Espoo, puh. 020 112 1510
asiakaspalvelu@mediq.com
www.mediq.fi/apuvaline

etac[®]
Creating Possibilities