




# Spare part list – EN

## Molift Mover 300

**molift**<sup>®</sup>  
by Etac

SPL06301 - 2026-04-20 - Rev. 6.0



[www.etac.com](http://www.etac.com)



Fig.1




Fig.3




Fig.2




Fig.	Pos.	Title	Item no.
1	1	Front wheel kit ø75 Mover 205/300/Partner 255	0950210
	2	Rear wheel ø100 w. brake Mover 205/300/Partner 255	0950211
	3	Locking handle kit chassis Mover 205/300/Partner 255	0950249
	4	Chassis contact kit Mover 205/300/Partner 255	0950246
	5	Special tool Molift Smart 150/Mover/Partner 255 (no illustration)	0920501
2	6	Leg adjust mechanism repair kit with motor	1420540
	7	Leg adjust mechanism repair kit without motor	1420530
	8	Motor, leg adjust mechanism	1320155
3	9	Actuator kit Mover 300	0950231
	10	Spare handle kit Mover 205/300/Partner 255	0950239
	11	Control box mounting kit	0950247
	12	Spare control box, BLE	0950252
	13	Channel1 cable f/BLE contr.box	1860116
	14	Channel2 cable f/BLE contr.box	1860117
	15	Spare actuator cable Mover 300	1860117
	16	Cable RJ50 f/handctrl. 2&4 btn (no illustration)	0950248

Fig.4




Fig.5




Fig.7




Fig.8




Fig.9



Fig.10



Fig.11



Fig.12




Fig.	Pos.	Title	Item no.
4	17	Cable for Molift Service Tool (USB/RJ50)	2510250
	18	Prog.box kit f/BLE control box (only for service technician)	0950254
5	19	Spare cable f/prog. box, RJ50 (only for service technician)	0950152
	20	USB A to USB B cable 1M (only for service technician)	0950153
6	21	Spiral cable USA plug (no illustration)	2920092
	22	Spiral cable EU plug (no illustration)	2920093
	23	Spiral cable GB plug (no illustration)	2920094
	24	Spiral cable AU plug (no illustration)	2920095
7		Product label	1860111
8		Label Safety Check	1100303
9		Label Safety Check	1110306
10		Label Molift by Etac	0950128
11		Label charging instructions (not available for Japan)	0950123
12		Label elec.data f/BLE ctrl.box (not available for Japan)	0950124



**Etac Ltd.**

Unit 60, Hartlebury Trading Estate,  
Hartlebury, Kidderminster,  
Worcestershire, DY10 4JB  
+44 121 561 2222



**Promeifa SA**

Avenue de Morges 33  
1027 Lonay, CH  
+41 44 872 97 79



**Etac Apm Media A/S**

Parallevej 3

DK-8751 Gedved



+45 796 858 33



patienthandling@etac.com



www.etac.com

**molift**<sup>®</sup>  
by Etac